



Recombinant Mouse Proto-oncogene Mas (Mas1)

Product Code	CSB-CF013505MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P30554
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MDQSNMTSLAEEKAMNTSSRNASLGSSHPPIPIVHWVIMSISPLGFVENGILLW FLCFRM RRNPFTVYITHLSIADISLLFCIFILSIDYALDYELSSGHHTIVTLVTFVFLGYNTG LY LLTAISVERCLSVLYPIWYRCHRPKHQSAFVCALLWALSCLVTTMEYVMCIDSG EESHRSR SDCRAVIIFIALSFLVFTPLMLVSSSTILVVKIRKNTWASHSSKLYIVIMVTIIIFLIFA MPMRVLYLLYYEYWSAFGNLHNISLLFSTINSSANPFIYFFVGSSKKKRFRESL KVVLTR AFKDEMQRPRQEGNGNTVSIETVV
Source	in vitro E.coli expression system
Target Names	Mas1
Protein Names	Recommended name: Proto-oncogene Mas
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.