



Recombinant Chicken Adenosine receptor A1 (ADORA1)

Product Code	CSB-CF001375CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49892
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MAQSVTAFQAAYISIEVLIALVSVPGNILVIWAVKMNQALRDATFCFIVSLAVADV AVGA LVIPLAIIINIGPQTEFYSCLMMPVLLITTESSILALLAIAVDRLRVKIPVRYKSV VT PRRAAVAIACCWIVSFLVGLTPMFGWNNLNKVLGTRDLNVSHSEFVIKCFET VISMAYM VYFNFFVWVLPPLLLMLLIYLEVFNLIRTQLNKKVSSSSNDPQKYYGKELKIAKS LALVL FLFALSWLPLHILNCITLFCPSCKTPHILTYIAIFLTHGNSAMNPVYAFRIKKFRT AFL QIWNQYFCCKTNKSSSSSTAETVN
Source	in vitro E.coli expression system
Target Names	ADORA1
Protein Names	Recommended name: Adenosine receptor A1
Expression Region	1-324
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.