



Recombinant Human Olfactory receptor 1S2 (OR1S2)

Product Code	CSB-CF844069HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGQ3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MKTLCSFLQISRNMHQENQTTITTEFILLGLSNQAEHQNLLFVLFVLSMYVVTVVG NGLIIV AISLDIYLHTPMYLFAYLSFADISSISNSVPKMLVNIQTNSQSISYESCITQMYFS IVF VVTDNLLLGTMAFDHFVAICHPLNYTTFMRRARFGTLLTVISWFLSNIIALHTLLLI QLL FCDHNTLPHFFCDLAPLLKLSCSDTMINELVLFIVGLSVIIFPFVLIFFSYVCIIRAV LG VSSTQGKWKAFSTCGSHLTIALLFYGTTVGVYFFPSSTHPEDTDKIGAVLFTV TPMMNP FIYSLRNKDMKGALRKLINRKISSL
Source	in vitro E.coli expression system
Target Names	OR1S2
Protein Names	Recommended name: Olfactory receptor 1S2 Alternative name(s): Olfactory receptor OR11-231
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.