



Recombinant Mouse Melanocortin receptor 5 (Mc5r)

Product Code	CSB-CF013562MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P41149
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MNSSSTLTVLNLTLNASEDGILGSNVKNKSLACEEMGIAVEVFLTLGLVSLLENI LVIGA IVKNKNLHSPMYFFVGLAVADMLVSMNAWETVTIYLLNNKHLVIADTFVRHI DNVFDS MICISVVASMCSLLAIAVDRIYITIFYALRYHHIMTARRSGVIIACIWTFCISCGIVFII Y YESKYVIIICLISMFFTMLFFMVSLYIHMFLARNHVKRIAASPRYNSVRQRTSMK GAILL TMLLGIFIVCWSPFFLHLILMISCPQNVYCSCFMSYFNMYLILIMCNSVIDPLIYAL RSQ EMRRTFKEIVCCHGFRRPCRLGGY
Source	in vitro E.coli expression system
Target Names	Mc5r
Protein Names	Recommended name: Melanocortin receptor 5 Short name= MC5-R
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.