



Recombinant Bovine Melanocortin receptor 5 (MC5R)

Product Code	CSB-CF013562BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P56451
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MNSSFHLHFLDLGLNTTDGNLSGLSVQNASSLCEDMGIAVEVFLALGLISLLENI LVIGA IVRNRNLHTPMYFFVGS LAVADMLVSLSNSWETITIIYLLTNKHLVMADASVRHL DNVFDS MICISVVASMCSLLAIAVD RYVTIFCALRYQRIMTGRRSGAIIGGIWAF CASC GT VFIVY YESTYVVICLIAMFLTMLLLMASLYTHMFL LARTHIRRIATLP GHSSVRQR TGVK GAILL AMLLGVFIVCWAPFFLHLILMISCPHNL YCSCFM SHFNMYLILIMCNSVIDPLIYA FRSQ EMRKT FKEIVCFQSFRTPCRFPSRY
Source	in vitro E.coli expression system
Target Names	MC5R
Protein Names	Recommended name: Melanocortin receptor 5 Short name= MC5-R
Expression Region	1-325
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.