



Recombinant Dog Gonadotropin-releasing hormone receptor (GNRHR)

Product Code	CSB-CF009637DO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9MZI6
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	MASASPEQNNHCSAVNNSNMLMQGNLPTLTLSGKIRVTVFFLFLSTIFNAS FLLKLQ KWTQKKEKGKLSRMKVLLKHLTLANLLETLIVMPLDGMWNITVQWYAGEFLC KVLSYLK LFSMYAPAFMMVVISLDRSLAITRPLAMKNNGKLGQSMIGLAWLLSGIFAGPQL YIFRMI HLADSSGQTEGFPQCVTHCSFPQWWHQAFFYNFFTFSCLFIIPLFITLICNAKIIFT LTRV LHQDPHELQLNQSKNNIPRARLRTLKMTVAFATSFTVCWTPYYVLGIWYWFDP EMLNRVS DPVNHHFFLFALLNPCFDPLIYGYSL
Source	in vitro E.coli expression system
Target Names	GNRHR
Protein Names	Recommended name: Gonadotropin-releasing hormone receptor Short name= GnRH receptor Short name= GnRH-R
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.