



Recombinant Human Olfactory receptor 9G4 (OR9G4)

Product Code	CSB-CF854140HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGQ1
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MIFPSHDSQAFTSVDMEVGNCTILTEFILLGFSADSQWQPILFGVFLMLYLITLS GNMTL VILIRTDSHLHTPMYFFIGNLSFLDFWYTSVYTPKILASCVSEDKRISLAGCGAQ LFFSC VVAYTECYLLAAMAYDRHAAICNPLLYSGTMSTALCTGLVAGSYIGGFLNAIAH TANTFR LHFCGKNIIDHFFCDAPPLVKMSCTNTRVYEKVLLGVVGFTVLSSILAILISYVNIL LAI LRIHSASGRHKAFSTCASHLISVMLFYGSLLFMYSRPSSTYSLERDKVAALFYT VINPLL NPLIYSLRNKDIKEAFRKATQTIQPQT
Source	in vitro E.coli expression system
Target Names	OR9G4
Protein Names	Recommended name: Olfactory receptor 9G4 Alternative name(s): Olfactory receptor OR11-216
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.