



Recombinant Human Olfactory receptor 6A2 (OR6A2)

Product Code	CSB-CF016993HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O95222
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MEWRNHSGRVSEFVLLGFPAPAPLQVLLFALLLLAYVLVTENTLIIMAIRNHST LHKPM YFFLANMSFLEIWYVTVTIPKMLAGFVGSKQDHGQLISFEGCMTQLYFFLGLGC TECVLL AVMAYDRYMAICYPLHYPVIVSGRLCVQMAAGSWAGGFGISMVKVFLISGLSY CGPNIIN HFFCDVSPLLNLSCTDMSTAELTDFILAIFILLGPLSVTGASYVAITGAVMHIPSA AGRY KAFSTCASHLTVVIIFYAASIFIYARPKALSAFDTNKLVSPLYAVIVPLLNPIIYCLR NQ EVKRALCCTLHLYQHQPDPKKASRV
Source	in vitro E.coli expression system
Target Names	OR6A2
Protein Names	Recommended name: Olfactory receptor 6A2 Alternative name(s): Olfactory receptor 11-55 Short name= OR11-55 Olfactory receptor 6A1 Olfactory receptor OR11-83 hP2 olfactory receptor
Expression Region	1-327
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.