



Recombinant Rat P2Y purinoceptor 6 (P2ry6)

Product Code	CSB-CF737185RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q63371
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MERDNGTIQAPGLPPTTCVYREDFKRLLLPPVYSVVLVVGLPLNVCVIAQICAS RRTLTR SAVYTLNLALADLLYACSLPLLIYNYARGDHWPFGDLACRLVRFLFYANLHGSIL FLTCL SFQRYLGICHPLAPWHKRGRRAAWVVCVWLVVTAQCLPTAVFAATGIQR NRTVCYDL SPPILSTRYLPYGMALTVIGFLLPFTALLACYCRMARRLCRQDGPAGPVAQER RSKAARM AVVVAAVFVISFLPFHITKTAYLAVRSTPGVSCPVLETFAAAYKGTRPFASANSV LDPIL FYFTQQKFRRQPHDLLQKLTAKWQRQRV
Source	in vitro E.coli expression system
Target Names	P2ry6
Protein Names	Recommended name: P2Y purinoceptor 6 Short name= P2Y6
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.