



Recombinant Human Olfactory receptor 4A16 (OR4A16)

Product Code	CSB-CF839911HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NH70
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MRPSSNVTEFVLLGLTQDPDVKKTLFVMFLLIYIVTMVGNLLIWVTTIGSPSLGS LMYFF LAYLSLMDAIYSTAMSPKLMIDLLCDKIAISLSACMGQLFIEHLLGGAEVFLLVVM AYDR YVAISKPLHYLNIMNRLVCILLLVAMIGGFVHSVQIVFLYSLPICGPNVIDHSV CDMY PLLELLCLDTYFIGLTVVANGGIICMVIFTFLLISCGVILNFLKTYSQEERHKALPT CIS HIIVVALVFVPCIFMYVRPVSNFPFDKLMTVFYSIITMLNPLIYSLRQSEMKNAM KNLW CEKLSIVRKRVSPTLNIFIPSSKATNRR
Source	in vitro E.coli expression system
Target Names	OR4A16
Protein Names	Recommended name: Olfactory receptor 4A16 Alternative name(s): Olfactory receptor OR11-117
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.