




Recombinant Human Olfactory receptor 13A1 (OR13A1)

Product Code	CSB-CF854141HU
Abbreviation	Recombinant Human OR13A1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q8NGR1
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose
Product Type	Transmembrane Proteins
Immunogen Species	Homo sapiens (Human)
Sequence	MKLWMESHLIVPETRPSRMMNSQTLVTEFILQGFSEHPEYRVFLFSCFLFLY SGALTGN VLITLAITFNPGLHAPMYFFLLNLATMDIICTSSIMPKALASLVSEESSISYGGCM AQLY FLTWAASSELLLLTVMAYDRYAAICHPLHYSSMMSKVFCGLATAVWLLCAVN TAHTGL MLRLDFCGPNVIIHFFCEVPPLLLLSCSSTYVNGVMIVLADAFYGIVNFLMTIASY GFIV SSILKVKTAWGRQKAFSTCSSHLTVVCMYYTAVFYAYISPVSGYSAGKSKLAG LLYTVLS PTLNPLIYTLRNKEVKAALRKLPFFRN
Source	in vitro E.coli expression system
Target Names	OR13A1
Protein Names	Recommended name: Olfactory receptor 13A1 Alternative name(s): Olfactory receptor OR10-3
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Protein Length	full length protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.