



Recombinant Rabbit Adenosine receptor A1 (ADORA1)

Product Code	CSB-CF001375RB
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P34970
Product Type	Transmembrane Protein
Immunogen Species	Oryctolagus cuniculus (Rabbit)
Sequence	MPPSISAFQAAYIGIEVLIALVSVPGNVLVIWAVKVNQALRDATFCFIVSLAVADV AVGA LVIPLAILINIGPETYFHTCLMVACPVLILTQSSILALLAIAVDRYLRVKIPLRYKAV VT PRRAAVAIAGCWILSLVVG LTPMFGWNNLREVQRAWAANGSVGEPVIKCEFE KVISM EYM VYFNFFVWVLP LLLMVLIIYLEV FYLIR RQLSKKASASSGDPHKYYGKELKIAKS LALIL FLFALSWLPLHILNCVTLFCPSCQKPSILVYTAIFLTHGNSAMNPIVYAFRIHKFR VTFL KIWN DHFR CRPAPAGDGD EDLPEEK PND
Source	in vitro E.coli expression system
Target Names	ADORA1
Protein Names	Recommended name: Adenosine receptor A1
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.