



Recombinant Pig Gonadotropin-releasing hormone receptor (GNRHR)

Product Code	CSB-CF009637PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P49922
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MANSASPEQNNQHCSAINSSILLTQGNLPTLTLSGKIRVTVFFLFLLSTAFNAS FLLKL QKWTQRKEKGGKLSRMKVLLKHLTLANLLETIVMPLDGMWNITVQWYAGEFL CKVLSYL KLFSMYAPAFMMVVISLDRSLAITRPLAVKSNSRLGRFMIGLAWLLSSIFAGPQL YIFRM IHLADSSGQTEGFSQCVTHGSFPQWWHQAFYNFFTFSCLFIIPLLIMLICNAKIM FTLTR VLQQDPHNQLNQSNNIPRARLRTLKMTVAFAASFIVCWTPYYVLGIWYWFD PEMVNRV SDPVNHFFFLFAFLNPCFDPLIYGYSFL
Source	in vitro E.coli expression system
Target Names	GNRHR
Protein Names	Recommended name: Gonadotropin-releasing hormone receptor Short name= GnRH receptor Short name= GnRH-R Alternative name(s): Luteinizing hormone-releasing hormone receptor Short name= LHRH
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.