



Recombinant Rat Free fatty acid receptor 2 (Ffar2)

Product Code	CSB-CF008605RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q76EI6
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MTPDWHSSLILTAYILIFLTGLPANLLALRAFVSRVRQPQPAPVHILLNLTLADL LLLL LLPFRIVEAASNFRWYLPKIVCALTGFGFYSSIYCSTWLLAGISIERYLGVAFPV QYKLS RRPLYGVIAALVAWIMSFHCTIVIIVQYLNSTEQVGTENQITCYENFTQAQLDV VLPVR LELCLVFFVPMTVTIFCYWRFVWIMLTQPHVGAQRRRRAVGLAVVTLNFLV CFGPYNM SHLVGFHLRQSPSWRVEAVVFSSLNASLDPLLFFYSSSVVRRAFGKGLLLLRN PGSSMLG RGAETVEGTKTDRGGSQTEGAQSSDFVTE
Source	in vitro E.coli expression system
Target Names	Ffar2
Protein Names	Recommended name: Free fatty acid receptor 2 Alternative name(s): G-protein coupled receptor 43
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.