



# Recombinant *Trachypithecus francoisi* Mas-related G-protein coupled receptor member X2 (MRGPRX2)

<b>Product Code</b>	CSB-CF680968TOZ
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q4QXU4
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	<i>Trachypithecus francoisi</i> (Francois' leaf monkey) ( <i>Presbytis francoisi</i> )
<b>Sequence</b>	MDPTTLVWGTESTTMNGNDQALPLLCGKETLILVVLILFIALVGLVGNFVLWLL GFRMR RNAFSVYVLSLAGADFLFLCFPMINCLAYLINFHISISINFPSFFTTVMTCAYL LAGLSML SAISTERCLSVLWPIWYRSRRPRHLSAVMCVLLWALSLLLSILEGKFCGFLFSD GDGWC QTFDFITAAWLMFLFVVLGSSLLALLVRILCGSRGLPLTRLYLTILLTVLIFLLCGL PFG IQWFLILWIWKNSDVLFCIHHPVSVLSSFNSSANPIIYFFVGSFRKQWRLRQP VKLALQRALQDTAEVDHSEGCFSQGTLEMSGSSLV
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	MRGPRX2
<b>Protein Names</b>	Recommended name: Mas-related G-protein coupled receptor member X2
<b>Expression Region</b>	1-330
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.