



Recombinant Pygathrix bieti Mas-related G-protein coupled receptor member X2 (MRGPRX2)

Product Code	CSB-CF683932FHN
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q4QXU2
Product Type	Transmembrane Protein
Immunogen Species	Rhinopithecus bieti (Black snub-nosed monkey) (Pygathrix bieti)
Sequence	MDPTTPAWGTESTTMNGNDQALPLLCGKETLILVLLILFIALVGLVGNFVLWLL GFRMR RNAFSVYVLSLAGADFLFLCFPMINCLEYLINFFHSISINFPSFFTTVMTCAYLAG LSML SAISTERCLSVLWPIWYRCRRPRHLSAVLCVLLWALSLLLSILEGKFCGLLFS DSDSGWC QTFDFITAAWLMFLFVVLGSSLLVLRILCGSQGLPLTRLYLTILLTVLIFLLCGL PFG IQWFLILWIWKNSDVLFCIHHPVSVLSSFNSSANPIIYFFVGSFRKQWRLRQP VLKAL QRALQDTAEVDHSEGCFSQGTLEMSGSSLV
Source	in vitro E.coli expression system
Target Names	MRGPRX2
Protein Names	Recommended name: Mas-related G-protein coupled receptor member X2
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.