



Recombinant *Macaca fascicularis* Mas-related G-protein coupled receptor member X2 (MRGPRX2)

Product Code	CSB-CF725943MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q5U9D9
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MDPTTPAWGTESTTMNGNDQALPLLCGKETMISVFLILFIALVGLVGNAFVLWL LGFRMR RNAFSVYVLSLAGADFLFLCFQMTNCLAYLINFFGSISINFPSFFTTVMTCACLA GLSML SAISTQRCLSILWPIWYRCRRPRHLSAVMCVLLWALSLLLSILEGKFCGFLFSD GDSGWC QTFDFITAAWLMFLFVVLGSSLALLVRILCGSRSLPLTRLYLTILLTVLIFLLCGL PFG IQWFLILWIWKNSDVLFCIHHPVSVVLSFNSSANPIIYFFVGSFRKQWRLRQPI LKLAL QRALQDTAEVDHSEGCFSQGTLEMSRSSLV
Source	in vitro E.coli expression system
Target Names	MRGPRX2
Protein Names	Recommended name: Mas-related G-protein coupled receptor member X2
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.