



Recombinant Bos indicus Duffy antigen/chemokine receptor (DARC)

Product Code	CSB-CF769670BUL
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q863H8
Product Type	Transmembrane Protein
Immunogen Species	Bos indicus (Zebu)
Sequence	MGNCLYPVADDNSTKLAIKEDLLIDFPEDYYPDYNEDVEAAAPCHSCSLLNYS SLPFFI LVSILGILASGTILYALLRPLFRWQLYQDRSTLVQLAVGSALFSIVVPILARGLSG ALIT SLCHLAHLVAYGSAFAQALLIGYHACLGPQLGAGQVPGLRLGVTVGLWGVAAL LSLPVVL GSDTSQGLCTVTFSGEWETLRYIHAAACFAIFVLLPLGLLGTKGLKTVLGRAPC PWVDVL WWFIFWWPQGMTLGLDSLVRSKAIVVSTCPAQQALDMLLDVAEALAILHCVA TPLLLAW VCYQATHTSPPSLPLPTTQTSHLDTLGSKS
Source	in vitro E.coli expression system
Target Names	ACKR1
Protein Names	Recommended name: Duffy antigen/chemokine receptor Alternative name(s): CD_antigen= CD234
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.