



Recombinant Human Olfactory receptor 2D3 (OR2D3)

Product Code	CSB-CF839890HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGH3
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MCSFFLCQTGKQAKISMGEENQTFVSKFIFLGLSQDLQTQILLFILFLIYLLTVLG NQL IIILIFLDSRLHTPMYFFLRNLSFADLCFSTSIVPQVLVHFLVKRKTISFYGCMTQII VF LLVGCTECALLAVMSYDRYVAVCKPLYYSTIMTQRVCLWLSFRSWASGALVSL VDSFTF HLPYWGQNIINHVFCEPPALLKLASIDTYSTEMAIFSMGVILLAPVSLILGSYW NIIST VIQMQSSEGRLKAFSTCGSHLIVVLFYSGGIFTYMRPNSKTTKELDKMISVFY TAVTPM LNPIIYSLRNKDVKGALRKLVGKCFSHRQ
Source	in vitro E.coli expression system
Target Names	OR2D3
Protein Names	Recommended name: Olfactory receptor 2D3 Alternative name(s): Olfactory receptor OR11-89
Expression Region	1-330
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.