



Recombinant Rat Proteinase-activated receptor 3 (F2r12)

Product Code	CSB-CF838728RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q920E1
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	SFNGNEHSFEEFPLSDIEGWTGATTTIKAKCPEESITTLHVNNATMGYLRSSLS TKVIPA IYILVFIGVPANIVTLWKLSSRTKSICLVIFHTNLAIADLLFCVTLFPFKIAYHLNGN DW VFGEVMCRVTTVAFYGNMYCAILILTCMGINRYLATVHPFTYRKLPKRNFTLLM CGVWWV MVVLYMLPLAILKQEYHLVQPGITTCHDVHDTCESPLPFQFYFVSLAFFGFLIP FVVSV FCYTTLIHKLNAQDRKWLRYIKAVLLILVIFTICFAPTNILIIHHANYYSNTDSL YFM YLIALCLGSLNSCLDPFLYFIMSKIVDQLTS
Source	in vitro E.coli expression system
Target Names	F2r12
Protein Names	Recommended name: Proteinase-activated receptor 3 Short name= PAR-3 Alternative name(s): Coagulation factor II receptor-like 2 Thrombin receptor-like 2
Expression Region	38-368
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.