



Recombinant Bovine Neuropeptides B/W receptor type 1 (NPBWR1)

Product Code	CSB-CF847579BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8MJV3
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MHNASYWGPERANTSCPAPAPTLGCPNASGPAPPLPPPLAVAVPVVYAVICA VGLAGNSA VLFVLLRAPRRKTVTNLFILNLAVADELFTLVPPVNIADFLRRWPFGELLCKLV VAVDQ YNTFSSLYFLTVMSADRYLVVLATAESRRVAGRTYGAARAVSLAVWGVATLVV LPFAVFA RLDEEQGRRQCVLVFPQPEALWWRASRLYTLVLGFAIPVSTICVLYTSLLCRLR AIRLDS HAKALDRAKKRVTVLVAILAVCLLVWTPYHLSTVVALTTDLPQTPLVIAVSYFIT SLSY ANSCLNPFLYAFLDDSFRRSLRQLLACRTTS
Source	in vitro E.coli expression system
Target Names	NPBWR1
Protein Names	Recommended name: Neuropeptides B/W receptor type 1 Alternative name(s): G-protein coupled receptor 7
Expression Region	1-331
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.