



Recombinant Serpentine receptor class gamma-3 (srg-3)

Product Code	CSB-CF344391CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P46572
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MSYEHGISGYTNFNENIHYFYQFAYLFTAICINRILYVIWVSQRHFYRNQSFYN LYSVD CFTSVLAMSNELIFTRSFLYFPQLCVSFSEIVKNSPVFMRIYYCLLSYLIAIKPVIH IFI AVNRMSCVMFPVTYSQNWSQKLRIMLIVIFLAPFLVIWNVLISDNFIGYVNGGF GISYTR RVTWASLSLMQFTLIILTTLITMVTMTTTFYKMTTMKKRIKASERALCIAAALISV GFL EAITQSFFAFFKEAPWLLDVMNYLRFATMDILFVGSPLVLLLVSQDQFRGHVLS RIGRTQ RVSSINNTSHIHHNTHHTMTRYSYFLWNVNK
Source	in vitro E.coli expression system
Target Names	srg-3
Protein Names	Recommended name: Serpentine receptor class gamma-3 Short name= Protein srg-3
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.