



Recombinant Rat Melanocortin receptor 4 (Mc4r)

Product Code	CSB-CF013561RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P70596
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MNSTHHHGMYSLHLWNRSSHGLHGNASESLGKGHSDGGCYEQLFVSPEVF VTLGVISLL ENILVIVAIKKNLHSPMYFFICSLAVADMLVSVSNGSETIVITLLNSTDTDAQS FTVN IDNVIDSVICSSLLASICSLLSIAVDRYFTIFYALQYHNIMTVRRVGIISCIWAACT VS GVLFIISDSSAVIICLITMFFTMLVLMASLYVHMFLMARLHIKRIAVLPGTGTRQ GAN MKGAITLTILIGVFVVCWAPFFLHLLFYISCPQNPYCVCFMSHFNLYLILIMCNAV IDPL IYALRSQELRKTFFKEIICFYPLGGICELPGRY
Source	in vitro E.coli expression system
Target Names	Mc4r
Protein Names	Recommended name: Melanocortin receptor 4 Short name= MC4-R
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.