



Recombinant Rat Adenosine receptor A2b (Adora2b)

Product Code	CSB-CF001377RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P29276
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	<p>MQLETQDALYVALELVIAALAVAGNVLVCAAVGASSALQTPTNYFLVSLATADV AVGLFA IPFAITISLGFCTDFHSCFLACFVLVLTQSSIFSL LAVAVDRYLAI RVP LRYKGLV TGT RARGIIAVLWVLAFGIGLTPFLGWNSKDRATSNCTEPGDGITNKSCCPVKCLFE NVVPM S YMVYFNFFGCVLPPLLIMMVIYIKIFMVACKQLQHMELM EHSRTTLQREIHA AK SLAMIV GIFALCWLPVHAINCITLFHPALAKDKPKWVMNVAILLSHANSVVNPIVYAYRNR DFRYS FHRIISRYVLCQTDTKGGSGQAGGQSTFSLSL</p>
Source	in vitro E.coli expression system
Target Names	Adora2b
Protein Names	Recommended name: Adenosine receptor A2b
Expression Region	1-332
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.