



Recombinant Pan paniscus Probable G-protein coupled receptor 33 (GPR33)

Product Code	CSB-CF670999EQU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q49SQ2
Product Type	Transmembrane Protein
Immunogen Species	Pan paniscus (Pygmy chimpanzee) (Bonobo)
Sequence	MDLINSTDYLINASTLVRNSTQFLAPASKMIIALSLYISSIIGTITNGLYLWVLRFK MKQ TVNTLLFFHLILSYFISTMILPFMATSQIQDNHWNFGTALCKVFNGTSLGLMFTS VFFLS AIGLDRYLLTLHPVWSQQHRTPRWASSIVLGVWISAAALSIPYLIFRQTHHDRK GKVTCQ NNYAVSTNWESEKEMQALRQWIHVACFISRFLGFLPPFFIIIFCYERVASKVKER SLFKS SKPFKVMMTAIIISFFVCWMPYHIHQGLLLTMNQSLLELTLILTTLTTSFNTIFSP TLYL FVGENFKKVFKKSILALFESTFSEDSSVERTQT
Source	in vitro E.coli expression system
Target Names	GPR33
Protein Names	Recommended name: Probable G-protein coupled receptor 33
Expression Region	1-333
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.