



Recombinant Serpentine receptor class alpha-13 (sra-13)

Product Code	CSB-CF630295CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q20618
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MAIISSVNRTCASESLELYRSYKYILSTSFNIIPIISLFFLVYAIKQLCAQSIIQYST RVLLITTILFAVCHQIAYFCFKADLLYTMLFKLDQPCNLQHSSYDCRFITIATTTT NCGM ALVQLAMSIDRVFALKFNRYVYKLSIPGITLALITLSISFSMFFILTIDDPLSGYVN HC GFYPTYSQDKFHIFLDVTLYLAVFNFDIGLMYYSYQEILWKRSYSFVNRFQS RISLKC TQAIFIISICQCISNVLYSGLLSLLMKLGRYMSSADYNLSLSLAYTTPYSCILPILI CK VLEYIKKQRTVGILSLRNQKQSMEGHMAMINSAWK
Source	in vitro E.coli expression system
Target Names	sra-13
Protein Names	Recommended name: Serpentine receptor class alpha-13 Short name= Protein sra-13
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.