



Recombinant Human Olfactory receptor 10R2 (OR10R2)

Product Code	CSB-CF818808HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8NGX6
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MPQILIFTYLNMFYFFPPLQILAENLTMVTEFLLLGFSSLGEIQLALFVVFLFLYLV ILS GNVTIISVIHLDKSLHTPMYFFLGILSTSETFYTFVILPKMLINLLSVARTISFNCCA LQ MFFFLGFAITNCLLLGVMGYDRYAAICHPLHYPTLMSWQVCGKLAACAIGGF LASLTVV NLVFSLPFCSANKVNHYFCDISAVILLACTNTDVNEFVIFICGVLVLVVPFLFICV SYLC ILRTLKIPSAEGRRKAFSTCASHLSVVIVHYGCASFIYLRPTANYVSNKDRLVTV TYTI VTPLLNPMVYSLRNKDVQLAIRKVLGKKGSLKLYN
Source	in vitro E.coli expression system
Target Names	OR10R2
Protein Names	Recommended name: Olfactory receptor 10R2 Alternative name(s): Olfactory receptor OR1-8
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.