



Recombinant Macaca mulatta Duffy antigen/chemokine receptor (DARC)

Product Code	CSB-CF853305MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95LF2
Product Type	Transmembrane Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Sequence	MGNCLHPAELSPSTQNSSQLNSDLWNFSYDGNDSFPDVDYDANLEAAAPCH SCNLLDDSA LPFFILVSVLGILASGTVLFMFFRPLFWQLCPGWPVLAQLAVGSALFSIVVPIL APGLG NTRSSALCSLGYCVWYGSAFAQALLGCHASLGPKLGAGQVPGLTLGLSVGL WGVAALLT LPITLASGASGGLCTPAYSMELKALQATHAVACLAVFVLLPLGLFGAKGLKKAL GMGPGP WMNILWAWFIFWWPHGVVGLDFLVRSKLLLLSTCLAQQALDLLLNLAELAIL HCVATP LLLALFCHQATRLLPSLPLPEGWSSHLDTLGSES
Source	in vitro E.coli expression system
Target Names	ACKR1
Protein Names	Recommended name: Duffy antigen/chemokine receptor Alternative name(s): CD_antigen= CD234
Expression Region	1-335
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.