



Recombinant Papio hamadryas Duffy antigen/chemokine receptor (DARC)

Product Code	CSB-CF822099ERJ
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95LG5
Product Type	Transmembrane Protein
Immunogen Species	Papio hamadryas (Hamadryas baboon)
Sequence	MGNCLHPAELSPSTQNSSQLNSEDLWNFSYDGNDSFPDVDYDANLEAAAPC HSCNLLDDS ALPFFILVSVLGILASGIVLFMFFRPLFWQLCPGWVLAQLAVGSALFSIVVPIL APGL GNTRSSALCSLGYCVWYGSFAQAALLGCHASLGPKLGADQVPGLTLGLSVG LWGVAALL TLPVTLASGASGGLCTPVYSMELKALQATHAVACLAIFVLLPLGLFGAKGLKKA LGMGPG PWMNILWAWFIFWWPHGVVLGLDFLVRSKLLLLSTCLAQQALDLLLNLAELAI LHCVAT PLLLALFCHQATRLLPSLPLPEGWSSHLDTLGSKS
Source	in vitro E.coli expression system
Target Names	ACKR1
Protein Names	Recommended name: Duffy antigen/chemokine receptor Alternative name(s): CD_antigen= CD234
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.