



Recombinant Rat 2-oxoglutarate receptor 1 (Oxgr1)

Product Code	CSB-CF754478RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q6Y1R5
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MIETLDSPANDSDFLDYITALENCTDEQISFKMQYLPVIYSIIFLVGFPNGNTVAISI YVF KMRPWKSSTIIMLNLALTDLLYLTSPLFLIHYYASGENWIFGDFMCKFIRFGFHF NLYSS ILFLTCFSLFRYIVIIHPMSCFSIQKTRWAVVACAGVWVISLVAVMPMTFLITSTT RTNR SACDLTSSDDLTTIKWYNLILTATTFCPLLLIVTLCYTTTIISTLTHGPRTHSCFKQ KAR RLTILLLLVFYVCFLPFHILRVIRIESRLLSISCSIESHIHEAYIVSRPLAALNTFGNL L LYVVVSNNFQQAFCSAVRCKAIGDLEQAKKDSCSNP
Source	in vitro E.coli expression system
Target Names	Oxgr1
Protein Names	Recommended name: 2-oxoglutarate receptor 1 Alternative name(s): Alpha-ketoglutarate receptor 1 G-protein coupled receptor 80 P2Y purinoceptor 15 Short name= P2Y15
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.