



Recombinant Mouse G-protein coupled receptor 26 (Gpr26)

Product Code	CSB-CF804891MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8BZA7
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MNSWDAGLAGLLVGTIGVSLLSNGLVLLCLLHSADIRRQAPALFTLNLTGNNLL CTVVNM PLTLAGVVAQRQPAGDRLCRLAAFLDTFLAANSMLSMAALSIDRWVAVVFPLS YRAKMRL RDAAFMVAYTWLHALTFPATALALSWLGFHQLYASCTLCSSRPDERLRFVAVFT SAFHALS FLLSFIVLCFTYLKVLKVARFHCKRIDVITMQTLVLLVDIHPVSRERCCLEEQKRRR QRAT KKISTFIGTFLVCFAPYVITRLVELFSTAPIGSHWGVLSKCLAYSKAASDPFVYSL LRHQ YRRSCKELLNRIFNRRSLHSVGLTGDSHSQNILPVSE
Source	in vitro E.coli expression system
Target Names	Gpr26
Protein Names	Recommended name: G-protein coupled receptor 26
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.