



Recombinant Human Visual pigment-like receptor peropsin (RRH)

Product Code	CSB-CF020517HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O14718
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLRNNLGNSSDSKNEGDSVFSQTEHNIVATYLMAGMISIISNIIVLGIFIKYKELR TPT NAIINLAVTDIGVSSIGYPMSAASDLYGSWKFGYAGCQVYAGLNIFFGMASIGL LTVVA VDRYLTICLPDVGRRMTTNTYIGLILGAWINGLFWALMPIIGWASYAPDPTGAT CTINWR KNDRSFVSYTMTVIAINFIVPLTVMFYCYHVTLSIKHHTTSDCTESLNRDWS QIDVTK MSVIMICMFLVAWSPYSIVCLWASFGDPKKIPPPMAIIAPLFAKSSTFYNPCIYV VANKK FRRAMLAMFKCQTHQTMPVTSILPMDVSQNPLASGRI
Source	in vitro E.coli expression system
Target Names	RRH
Protein Names	Recommended name: Visual pigment-like receptor peropsin
Expression Region	1-337
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.