



Recombinant Serpentine receptor class delta-32 (srd-32)

Product Code	CSB-CF522320CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O01608
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MFVDIYDQLYISYAAVVSTLGIIFNGFLLFLIFFKSPSCLTPYTVFLANTSITQLGY CIC FLLTVPRVISINLRIVNIYLGFSQFLGHWWSYMIFTTMLHFAVNSFLSIMLSMVFR CISL KTLRFPTSGAFAMCILAYMIPLSMVVSIRGIEITSNFTINSKYTLWQLENLDKYRT VVG TMAQLSTLWVACCVSILCIPYISVMFWCRYRILRMLERPGYMFNTTTTLQIKRL VKALTV QSLIPVFTLFPASLIFLSTQFHVIETTKFGYIIISLLSLSPTIDPLVTIYYVQPYRKYI V DLLWSEERPMVSPFLSNNDPRFYRSRSNSVLMRNTHFV
Source	in vitro E.coli expression system
Target Names	srd-32
Protein Names	Recommended name: Serpentine receptor class delta-32 Short name= Protein srd-32
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.