



Recombinant Rat Cysteinyl leukotriene receptor 1 (Cysltr1)

Product Code	CSB-CF006465RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q924T8
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MVGAENMTASFNNRCHDTIDEFRNQVYSTMYSMISVVGFFGNSFVLYVLIKT YHEKSAF QVYMINLAIADLLCVCTLPLRVVYYVHKGKWFFGDFLCRLTTYALYVNLYCSIFF MTAMS FFRCVAIVFPVQNINLVTQKKARFVCVGIWIFVILTSSPFLLSKSYQDEKNNTKC FEPPQ DKQTKKYVLVLHYVSLIFGFIPFVTIIVCYTMIILTLTKNTMCKNLPSRRKAIGMIIV V TAAFLVSFMPYHIQRAIHLHLHSETRSCDSVLRMQKSVVITLSLAASNCCFDP LLYFFS GGNFRRLSTFRKHSLSSMTYIPKKKASLPEKGEEMCKE
Source	in vitro E.coli expression system
Target Names	Cysltr1
Protein Names	Recommended name: Cysteinyl leukotriene receptor 1 Short name= CysLTR1
Expression Region	1-339
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.