



Recombinant Archaeoglobus fulgidus UPF0118 membrane protein AF_1800 (AF_1800)

Product Code	CSB-CF518054DOC
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O28474
Product Type	Transmembrane Protein
Immunogen Species	Archaeoglobus fulgidus (strain ATCC 49558 / VC-16 / DSM 4304 / JCM 9628 / NBRC 100126)
Sequence	MNKYWTLVLLLSILVVLALTFYFFTPLLDGIVMGVVFAYVAKPVKRKIEHVGRIK ASVIA TAIVILPISVLMFYGLIQGLNQAIYLITHYKVIESGILNILSKMGIEEGEEYVKWITS NV FSILQSSIQPSAVEITKKATLLILNFFISIVVCFYALADMENFVRRRTTSVVPEDRRD EFR RFVEEIDVTFESLWFGNFVAILIGLVSLPFFLYFNVPFAPLLSGLMFLAALIPIFA EWM IILPVSLYLLLVDVGRGLSFLIGVVFLYVLPILRPYFVGYTSKIHPLVLMLAFIG GG LVGGISGFFIAPMIVGLATAIYNYTKEELATENHDPEVS
Source	in vitro E.coli expression system
Target Names	AF_1800
Protein Names	Recommended name: UPF0118 membrane protein AF_1800
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.