



Recombinant Pig Cysteinyl leukotriene receptor 1 (CYSLTR1)

Product Code	CSB-CF006465PI
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q95N02
Product Type	Transmembrane Protein
Immunogen Species	Sus scrofa (Pig)
Sequence	MDGVRNLTVSCASSNTCNDDTIDDFRNQVYSTLYSMITVVGFFGNGFVLYVLIKT YHEKSA YQVYMINLAVADLLCVCTPLRVVYYVHKGWLFGLCRLSTYALYVNLVYCSIF FMTAM SFFRCIAIVFPVQINLITHKKAKIVCIAIWIFVILTSSPFLMSTSYKDEKNNTKCFE PP QXNQAKYHVLVLHYVSLFVGFIIIPFVIIVCYTMILTLLKNSMKKNISSRKKKAIGMII V VTAAFLISFMPYHIQRTIHLHLHNDTKHCDSVLRMQKSVXITLSLAASNCCFDP LLYFF SGGNFREGLSTFRKHSLSTMTYVPKKKTSLPEKAQEIIYKE
Source	in vitro E.coli expression system
Target Names	CYSLTR1
Protein Names	Recommended name: Cysteinyl leukotriene receptor 1 Short name= CysLTR1
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.