



Recombinant *Macaca fascicularis* Melanin-concentrating hormone receptor 2 (MCHR2)

Product Code	CSB-CF013585MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8SQ54
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MNPFHSSCWNTSAELSNKSWNKEFAYQTASAVDTVILPSMIGIICSTGLVGNILIVFTII RSRKKTVPDIYICNLAVADLVHIIGMPFLIHQWARGGEWVFGGPLCTIITSLDTCNQFAC SAIMTVMSVDRYFALVQPFRLTSWRTRYKTIRINLGLWAASFILALPVWIYSKVIKFKDG VESCAFDLTSPDDVLWYTLTYLTITTTFFFPLPLILVCYILILCYTWEMYQQNKDARCCNPS VPKQRMKLTKMVLVAVFILSAAPYHVIQLVNLQMEQPTLAFYVGYYSICLSYASS INPFLYILLSGNFQKRLPQIQRRVTDKEIKNMGNTLKSHF
Source	in vitro E.coli expression system
Target Names	MCHR2
Protein Names	Recommended name: Melanin-concentrating hormone receptor 2 Short name= MCH receptor 2 Short name= MCH-R2 Short name= MCHR-2 Alternative name(s): G-protein coupled receptor 145 GPRv17 MCH-2R Short name= MCH2
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.