



Recombinant Guinea pig Cysteinyl leukotriene receptor 1 (CYSLTR1)

Product Code	CSB-CF006465GU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q2NNR5
Product Type	Transmembrane Protein
Immunogen Species	Cavia porcellus (Guinea pig)
Sequence	<p>MDETGNPTIPPASNNTCYDSIDDFRNQVYSTLYSMISVVGFFGNGFVLYVLVKT YHEKSA FQVYMINLAVADLLCVCTPLRVAYYYVHKGIWLFGDFLCRLSTYALYVNLYCSIF FMTAM SFFRCVAIVFPVQNISLVTQKKARLVCAIWMFVILTSSPFLMANTYKDEKNNTK CFEPP QDNQAKNYVLILHYVSLFIGFIIPFITIIVCYTMIIFTLLKSSMKKNLSSRKRAIGMII V VTAAFLVSFMPYHIQRTIHLHFLHNKTKPCDSILRMQKSVVITLSLAASNCCFDP LLYFF SGGNFRRLSTIRKYSLSMITYIPKKKTSLPQKGKDICKE</p>
Source	in vitro E.coli expression system
Target Names	CYSLTR1
Protein Names	<p>Recommended name: Cysteinyl leukotriene receptor 1 Short name= CysLTR1 Alternative name(s): Cysteinyl leukotriene D4 receptor Short name= LTD4 receptor</p>
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>