



Recombinant Bovine Organic solute transporter subunit alpha (OSTA)

Product Code	CSB-CF664892BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q3T124
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MEPDRTQIRLDPRYTADLLEILKTNYSVPSACFSYPPTAAQLLRALGPVDISLMV IMTLF VLGSIAIFLEAAVYLHKNTRCPIKRKTLIWCSSSPTIVSAFSCFGLWIPRALTLVE MAIT TFYSMCFYLLMQAMVEGFGGKEAVLRTLKDPVMIHTGPCCCCCPCCPRIKIT RKRLQLL LLGPIQYAFFKISLTLVGLFLIPDGIFDPSDISEGSTALWINTFLGVSTLSALWTIGI IF RQARLHLGEGNIGAKFVLFQALLILSALQPSIFSVLASGGQIACSPPFSSKIRSQ VMNCH LLILESFLITVLTRIIYRRKDDKLGYPFSSPDQDLNLKA
Source	in vitro E.coli expression system
Target Names	SLC51A
Protein Names	Recommended name: Organic solute transporter subunit alpha Short name= OST-alpha
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.