



Recombinant *Macaca mulatta* (Rhesus macaque) C5a anaphylatoxin chemotactic receptor (C5AR1)

Product Code	CSB-CF003996MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79188
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca mulatta</i> (Rhesus macaque)
Sequence	TPDYGHYDDKDTLDANTPVDKTSNTLRVPDILALVIFAVVFLVGVLARNALVVVW TAFEAK RTINAIWFLNLAVADFLSCLALPILFTSIVQHHHWPFGGAACRILPSLILLNMYAS LLL ATISADRFLLVFNPIWCQNF RGAGLAWIACAVAWGLALLLTIPSFLYRVVREEY FPPKVL CGVDHGHDKRRERAVAIARLVLGFVWPLLTLTMCYTFLLLR TWSRRATRSTKT LKVVVAV VASFFIFWLPYQVTGMMMSFLEPSSPTFLLKLDLCLISFAYINCCINPIIYVVA GQGF QGRLRKSLPSLLRNVLTEESMVRESKSFTRSTVDTMAQKT
Source	in vitro E.coli expression system
Target Names	C5AR1
Protein Names	Recommended name: C5a anaphylatoxin chemotactic receptor Short name= C5a-R Short name= C5aR Alternative name(s): CD_antigen= CD88
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.