



Recombinant Chicken Adenosine receptor A2b (ADORA2B)

Product Code	CSB-CF001377CH
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O13076
Product Type	Transmembrane Protein
Immunogen Species	Gallus gallus (Chicken)
Sequence	MNTMKTTYIVLELIIAVLSIAGNVLCWAVAINSTLKNATNYFLVSLAVADIAVGLL AIP FAITISIGFQVDFHSCLEFFACFVLVLTQSSIFSLLAVIDRYLAIKIPLRYNSLVTGK RA RGLIAVLWLLSFVIGLTPLMGWNKAMSGCPNSTNETGADHGAGHHGCFISCLF ENVVTMS YMVYFNFFGCVLLPLIIMLGIYIKIFMVACKQLHQIELMGNSRTTLQKEVHAAKSL AIIV GLFAFCWLPLHILNCITHFHEEFSSKSKPEWVMYVAILSHANSVINPIIYAYRIRDF RYT FHKIISKILCKTDDFPKCTTDNNQHLLVTNVNAPAASVTI
Source	in vitro E.coli expression system
Target Names	ADORA2B
Protein Names	Recommended name: Adenosine receptor A2b Alternative name(s): Adenosine receptor 2b
Expression Region	1-340
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.