



Recombinant Goat Platelet-activating factor receptor (PTAFR)

Product Code	CSB-CF018946GO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9GK76
Product Type	Transmembrane Protein
Immunogen Species	Capra hircus (Goat)
Sequence	MEPNNSFRVDSEFRYTLFPIFYISIVFVLGVIANSYVLWVVFARLYPSKKFNEIKIF MVNLT MADLLFLVTLPLWIVYYYNQGDWILPKFLCNLAGCFFFINTYCSVAFLAVITYNR FQAVT RPIKTAQATTRKRGFLLSLIIVWSIVGAASYFFVLDSTNSEPKKTGSGNITRCFE HYEKG SIPVLIHIFLVFSFFLVFLIILFCNLVIIRTLTLLTQQVQMQRNAEVKRRALWMVCTVL AV FVICFVPHHLVQLPWTLAELGFQDQDFHQQGINDAHQVTLCLLSTNCVLDPIIYCF LTKKF RKHLTEKLYSMRESRKCSRATSETGTEVVVQLKDAPIKSLKY
Source	in vitro E.coli expression system
Target Names	PTAFR
Protein Names	Recommended name: Platelet-activating factor receptor Short name= PAF-R Short name= PAFr
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.