



Recombinant Bovine Platelet-activating factor receptor (PTAFR)

Product Code	CSB-CF018946BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9TTY5
Product Type	Transmembrane Protein
Immunogen Species	Bos taurus (Bovine)
Sequence	MEPNNSFRVDSEFRYTLFPIFYISIVFVLGVIANSYVLWVWFARLYPSKKFNEIKIF MVNLT MADLLFLVTLPLWIVYYYNQGDWILPKFLCNLAGCFFFINTYCSVAFLAVITYNR FQAVT RPIKTAQATTRKRGILLSLIWVSIVGAASYFFVLDSTNREPNTGSANITRCFEH YEKG SIPVLTIHIFLVFSFFLVFLIILFCNLVIIRTLTLLTQQVQIQRNAEVKRRALWMVCTVL AV FIICFVPHHLVQLPWTLAELGFQDQDFHQAINDAHQVTLCLLSTNCVLDPIIYCFL TKKF RKHLTEKLYSMRESRKCSRATSETGTEVVMQLKDVPVKSLKY
Source	in vitro E.coli expression system
Target Names	PTAFR
Protein Names	Recommended name: Platelet-activating factor receptor Short name= PAF-R Short name= PAFr
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.