



Recombinant Serpentine receptor class gamma-69 (srg-69)

Product Code	CSB-CF621175CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q19258
Product Type	Transmembrane Protein
Immunogen Species	Caenorhabditis elegans
Sequence	MNSCHPPQDEMAGLIGIYAFQGFYGLLSVVVYTFNIRALRHHKNNLDKSFSLLY TCCAAL SLTYFLDHFILIRRFVKLGFFCEIILENFGEPNYWMMPYKTIASYCPIAILVFHALIA AHR FSIVAAPMRGVQLWDRYRRLFVLVGFILPLIFMWFMI PCKSYAELDSEGSGGLD IEYKKV FSISSSLAAAIAAVLFGVLTCLCTFGMLIALAKLSLRKLSQAEISLIVFEVFM TVFT LIY AFTQGILYYSIYIVKDMELKSTVIQFRTFAIDIFILPQAWTLLFLSTTVRRYTLRAF GKR LGVEFLSTEIEKSARMVSVAPATISLQKSTVLNYNFTLQNLF
Source	in vitro E.coli expression system
Target Names	srg-69
Protein Names	Recommended name: Serpentine receptor class gamma-69 Short name= Protein srg-69
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.