



Recombinant Human Zinc transporter ZIP11 (SLC39A11)

Product Code	CSB-CF854023HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8N1S5
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MLQGHSSVFQALLGTTFTWGMTAAGAALVFVFSSGQRRILDGSLGFAAGVML AASYWSLL APAVEMATSSGGFGAFAFFPVAVGFTLGAAFVYLADLLMPHLGAAEDPQTTLA LNFGSTL MKKKSDPEGPALLFPESELSIRIGRAGLLSDKSENGEAYQRKKAAATGLPEGP AVPVPSR GNLAQPGGSSWRRIALLILAITIHNVPEGLAVGVGFGAIEKTASATFESARNLAI GIGIQ NFPEGLAVSLPLRGAGFSTWRAFWYGQLSGMVEPLAGVFGAFVLAEPILP YALAF AAG AMVYVVMDDIPEAQISGNGKLASWASILGFVVMMSLDVGLG
Source	in vitro E.coli expression system
Target Names	SLC39A11
Protein Names	Recommended name: Zinc transporter ZIP11 Alternative name(s): Solute carrier family 39 member 11 Zrt- and Irt-like protein 11 Short name= ZIP-11
Expression Region	1-342
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.