



Recombinant Mouse Formyl peptide receptor-related sequence 3 (Fpr-rs3)

Product Code	CSB-CF527121MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	O88537
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MEANSSIPLNGSEVVFYDSTTSRVLWILSVIVLSITFVLGVLGNGLVIWVAGFRM AHTVT TICYLNLALGDFSFMTLPLHIISMVMKGGKWLFGWFLCKFVLSIVHINLFSVFLI TLIA MDRCTCVLHPVWVQNHRTVSLARKVIVGAWILSLLLTLPHFLFTTVRDARGE VHCTCNF ESVVANPEEQQLKVSITVSTATGIISFIIGFSLPMSFIAVCYGLMAAKICRKGFLNS SRPL RVLTAVAISFFMCWFPPQLIILLGNIWNKETPSSIHILLNPASTLASFNSCLNPILY VFL GQEFREKLIYSLASLERALREDSVLSSGKSSNFSSCPADSEL
Source	in vitro E.coli expression system
Target Names	Fpr-rs3
Protein Names	Recommended name: Formyl peptide receptor-related sequence 3
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.