



Recombinant Mouse Mas-related G-protein coupled receptor member F (Mrgprf)

Product Code	CSB-CF844992MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q8VCJ6
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MAGNCSWEAHSTNQNKMCPLMSEARELYSRGFLTIEQIATLPPPAVTNYIFLL LCLCGLV GNGLVLWFFGFSIKRTPFSIYFLHLASADGMYLFSKAVIALLNMGTFLGSFPDYI RRVSR IVGLCTFFTGVSLIPAISIERCVSVIFPTWYWRRRPKRLSAGVCALLWMLSFLVT SIHNY FCMFLGHEAPGTVCRNMDIALGILLFFLCPLMVLPLCLALILHVECRARRRQRS AKLNHV VLAIVSVFLVSSIYLGIDWFLFWVFQIPAPFPEYVTDLCICINSSAKPIVYFLAGRD KSQ RLWEPLRVVFQRALRDGAEPGDAASSTPNTVTMEMQCPSGNAS
Source	in vitro E.coli expression system
Target Names	Mrgprf
Protein Names	Recommended name: Mas-related G-protein coupled receptor member F Short name= Mas-related gene F protein
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.