



Recombinant *Macaca fascicularis* C-X-C chemokine receptor type 6 (CXCR6)

Product Code	CSB-CF861038MOV
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9BDS6
Product Type	Transmembrane Protein
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Sequence	MAEYDHYEDDGFLNSFNDSSQEEHQDFLQFRKVFLPCMYLVVFCGLVGNL VLVISIFY HKLQSLTDVFLVNLPLADLVFVCTLPFWTYAGIHEWIFGQVMCKTLLGVYTINF YTSMLI LTCITVDRFIVVVKATKAYNQQAKRMTWGKVICLLIWWISLLVSLPQIIYGNVFNL DKLI CGYHDEEISTVVLATQMTLGFFLPLLAMIVCYSVIKTLHAGGFQKHRSLKIIFL VMAV FLLTQTPFNLVKLIRSTRWEYYAMTSFHYTIIVTEAIAYLRACLNPVLYAFVSLKF RKNF WKLVKDIGCLPYLGVSHQWKSSDNSKTFSASHNVEATSMFQL
Source	in vitro <i>E.coli</i> expression system
Target Names	CXCR6
Protein Names	Recommended name: C-X-C chemokine receptor type 6 Short name= CXC-R6 Short name= CXCR-6 Alternative name(s): G-protein coupled receptor STRL33 G-protein coupled receptor bonzo CD_antigen= CD186
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.