



Recombinant Human Thromboxane A2 receptor (TBXA2R)

Product Code	CSB-CF023260HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P21731
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	MWPNGSSSLGPCFRPTNITLEERRLIASPWFAASFCVVGLASNLLALSVLAGAR QGGSHTR SSFLTFLCGLVLTDFLGLLVTGTIVVSQHAALFEWHAVDPGCRLCRFMGVVMIF FGLSPL LLGAAMASERYLGITRPFSRPAVASQRRRAWATVGLVWAAALALGLLPLLGVGR YTVQYPG SWCFLTGAESGDVAFGLLFSMLGGLSVGLSFLNNTVSVATLCHVYHGQEAA QQRPRDSE VEMMAQLLGIMVVASVCWLPLLVFIAQTVLRNPPAMSPAGQLSRTTEKELLIYL RVATWN QILDPWVYILFRRRAVLRRQLPRLSTRPRSLSLQPQLTQRSGLQ
Source	in vitro E.coli expression system
Target Names	TBXA2R
Protein Names	Recommended name: Thromboxane A2 receptor Short name= TXA2-R Alternative name(s): Prostanoid TP receptor
Expression Region	1-343
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.