



Recombinant Mouse G-protein coupled receptor family C group 5 member D (Gprc5d)

Product Code	CSB-CF878065MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q9JIL6
Product Type	Transmembrane Protein
Immunogen Species	Mus musculus (Mouse)
Sequence	MYEDCVKSTEDYYLFCDNEGPWAIVLES LAVIGIVVTILLLLAFLFLMRKVQDCS QWNVL PTQFLFLLAVLGLFGLTFAFIIQLNHQTAPVRYFLFGVLF AICFSCLLAHASNLVK LVRG RVSF CWTTILFIAIGV SLLQTIIAIEYVTLIMTRGLMFEHMT PYQLNVDFVCLLIYV LFL MALTFFVSKATFCGPCENWKQHGR LIFATVLSIIIWVWVWISMLLRGNPQLQRQ PHWDDA VICIGLVTNAWVFLLIYIPELSILYRSCRQECPTQGNVCQVPVYQRSFRMDTQE PTRAR DSDGAQEDVALTAYGTPIQLQSADPSREYLIP SATLSPQQDAGL
Source	in vitro E.coli expression system
Target Names	Gprc5d
Protein Names	Recommended name: G-protein coupled receptor family C group 5 member D
Expression Region	1-344
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.